<p>Gebruiksaanwijzing</p> <p>OPHTEIS FR <sup>PRO</sup></p>	<p>Instruções de utilização</p> <p>OPHTEIS FR <sup>PRO</sup></p>	<p>Инструкция по применению</p> <p>OPHTEIS FR <sup>PRO</sup></p>	<p>Bruksanvisning</p> <p>OPHTEIS FR <sup>PRO</sup></p>	<p>Kullanım Talimatları</p> <p>OPHTEIS FR <sup>PRO</sup></p>
<p>NL</p>	<p>PT</p>	<p>RU</p>	<p>SV</p>	<p>TR</p>

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	<p>Manufacturer</p>
	<p>Use-by date</p>
	<p>Batch code</p>
	<p>Sterilized using ethylene-oxide (packaging and cannula)</p>
	<p>Sterilized using steam (product)</p>
	<p>Do not re-use</p>
	<p>Do not resterilize</p>
	<p>Caution</p>
	<p>Consult instructions for use</p>
	<p>Do not use if package is damaged</p>
	<p>Storage temperature 2<span> </span>°C - 25<span> </span>°C</p>
	<p>Keep dry</p>
	<p>Keep away from sunlight</p>
	<p>Does not contain latex</p>

	<p>For cannula</p> <p> Sterimedix Limited 1 Madeley Road, Redditch, B98 9NB, UK</p>
	<p></p>
	
<b>Customer Feedback</b>	
Please report any feedback, adverse event/ complaints to: <b>Rayner</b>	
The Ridley Innovation Centre 10 Dominion Way, Worthing West Sussex, BN14 8AQ, United Kingdom Tel.: +44 (0) 1903258900 Fax.: +44 (0) 1903258901 Email: feedback@rayner.com	

Ophteis FR Pro IFL PF000103 Version 02, Date of Issue: 1st November 2017		<b>The Ridley Innovation Centre</b> 10 Dominion Way, Worthing, West Sussex, BN14 8AQ, United Kingdom Tel: +44 (0)1903 258900 Web: rayner.com
--	---	--

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Instructions for use
OPHTEIS FR <sup>PRO</sup>

## EN

**Content:** Ophteis FR Pro is a viscoelastic sodium hyaluronate of non-animal origin and sorbitol solution. Each product consists of 1.2 mL of viscoelastic solution in a single-use glass syringe for intraocular use, together with one sterile 27 G cannula.

**Description:** Ophteis FR Pro is a viscoelastic, sterile, apyrogenic, isotonic solution of highly purified sodium hyaluronate of non-animal origin and sorbitol, in a physiological phosphate buffer. Sodium hyaluronate is a biopolymer obtained from bacterial biofermentation that belongs to the class of glycosaminoglycans. It possesses favorable physical and rheological properties for ophthalmic surgery. Sorbitol is well known as a safe excipient and it is a compound which is found in the human aqueous humor. 1 mL of Ophteis FR Pro contains 20 mg sodium hyaluronate, 40 mg sorbitol, sodium chloride, disodium hydrogen phosphate, disodium dihydrogen phosphate and water for injection. Ophteis FR Pro has a pH of 6.8 to 7.4 and osmolality of 295 to 355 mOsm/kg. The molecular weight of Ophteis FR Pro is 1.8 MDa. Ophteis FR Pro has a zero shear viscosity of 217-500 Pa s. The specific combination of high molecular weight sodium hyaluronate and sorbitol in Ophteis FR Pro allows an efficient lubrication, support and protection of ocular tissues during ophthalmic surgery. In particular, the free radical scavenging capacity obtained thanks to this combination offers a high protection during all the surgery of the intraocular tissues, such as the protection of the corneal endothelium. Due to its solubility in water and its rheological properties, Ophteis FR Pro can be easily removed of the intraocular space by irrigation and aspiration at the end of the surgery.

**Indications:** Ophteis FR Pro is an injectable solution for use as a surgical aid in surgery of the anterior and posterior segments of the eye, and especially:

- Cataract surgery with or without intraocular lens implantation
- Glaucoma surgery

**Contraindications:** Care should be taken in patients with hypersensitivity to any components in this solution.

**Directions for Use:** Ophteis FR Pro is carefully injected into the anterior chamber using the supplied sterile single-use cannula. It is recommended for intraocular lens implantation to cover the implant and instruments with Ophteis FR Pro immediately before surgery. This measure is to additionally protect the endothelium and the surrounding tissue. The Ophteis FR Pro volume to be injected varies and depends on the type of surgery. As a substitute for viscoelastic loss by flow or irrigation, Ophteis FR Pro can be injected several times. At the end of surgery, Ophteis FR Pro must be removed completely using a suitable irrigation/aspiration device.

**Warnings and Precautions:** Residual viscoelastic not removed at the end of the surgical procedure is flushed out naturally via the trabecular system and Schlemm’s canal. There is, however, a risk of blockage of the drainage channels, which can lead to an increase in intraocular pressure. To prevent the risk of IOP irrigation/aspiration is used at the end of the surgical procedure to remove residual viscoelastic. In case of IOP increase, treat with the preferred therapy. Patients who have been diagnosed with wide-angle glaucoma, severe myopia, diabetic retinopathy or uveitis are at particular risk of increased intraocular therapy.

- Do not reuse syringe and cannula. Any repeat usage of the syringe and cannula carries a risk of contamination and infection of the patient.
- Do not re-sterilize the pre-filled syringe and cannula.
- Do not use if package is damaged or opened.
- Do not use after the expiry date printed on the pack.
- Dispose of the syringe and cannula in accordance with accepted medical practice and applicable national, local and institutional requirements. There is no evidence concerning the safety of Ophteis FR Pro in human pregnancy and lactation. Administration during pregnancy and lactation is at the discretion of the ophthalmic surgeon.

**Storage:** Store between 2 °C – 25 °C. Protect from light and shocks. Do not freeze.

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