<b>Gebruiksaanwijzing</b> Methylvisc 2,0%	<b>NL</b>	<b>Instruções de utilização</b> Methylvisc 2,0%	<b>PT</b>	<b>Инструкция по применению</b> Methylvisc 2,0%	<b>RU</b>	<b>Bruksanvisning</b> Methylvisc 2,0%	<b>SV</b>	<b>Kullanım Talimatları</b> Methylvisc 2,0%	<b>TR</b>
--	-----------	--	-----------	--	-----------	--	-----------	--	-----------

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
	<b>1. Twist and pull out</b>
	<b>2. Insert and twist the cannula</b>
	<b>WRONG</b> (Backstop opening at front)
	<b>CORRECT!</b> (Backstop opening at back)

	<b>Manufacturer</b>
	<b>Use-by date</b>
	<b>Batch code</b>
	<b>Sterilized using ethylene-oxide (packaging and cannula)</b>
	<b>Sterilized using steam (product)</b>
	<b>Do not re-use</b>
	<b>Do not resterilize</b>
	<b>Caution</b>
	<b>Consult instructions for use</b>
	<b>Do not use if package is damaged</b>
	<b>Storage temperature 2<span> </span>°C - 25<span> </span>°C</b>
	<b>Keep dry</b>
	<b>Keep away from sunlight</b>
	<b>Does not contain latex</b>

	<b>For cannula</b> Sterimedix Limited 1 Madeley Road, Redditch, B98 9NB, UK

**Customer Feedback**  
Please report any feedback, adverse event/ complaints to:
**Rayner**  
The Ridley Innovation Centre  
10 Dominion Way, Worthing  
West Sussex, BN14 8AQ, United Kingdom  
Tel.: +44 (0) 1903258900  
Fax.: +44 (0) 1903258901  
Email: feedback@rayner.com

Methylvisc 2.0% IFU  
PF000105  
Version 02, Date of issue: 1st November 2017

# METHYLVISC 2.0

<b>EN</b>	Instructions for use
<b>AR</b>	تعليمات الاستخدام
<b>CS</b>	Pokyny pro použití
<b>DE</b>	Gebrauchsanweisung
<b>ES</b>	Instrucciones de uso
<b>FR</b>	Mode d'emploi
<b>HU</b>	Használati utasítás
<b>ID</b>	Petunjuk penggunaan
<b>IT</b>	Istruzioni per l'uso
<b>KO</b>	사용 지침서
<b>NL</b>	Gebruiksaanstructies
<b>PT</b>	Instruções de utilização
<b>RU</b>	Инструкция по применению
<b>SV</b>	Bruksanvisning
<b>TR</b>	Kullanım talimatları

**The Ridley Innovation Centre**  
10 Dominion Way, Worthing,  
West Sussex, BN14 8AQ, United Kingdom  
Tel: +44 (0)1903 258900  
Web: rayner.com

Instrukcije za upotrebu	EN	ES	FR	HR	IT	KO
Methylvisc 2.0%	تعليمات الاستخدام	Metilvisc 2,0%	Metilvisc 2,0%	Metilvisc 2,0%	Metilvisc 2,0%	메틸비스크(Methylvisc) 2.0%

**Content:** Methylvisc 2.0% is a viscoelastic Hydroxypropyl Methylcellulose solution.

Each product consists of 2.0 mL of viscoelastic solution in a single-use glass syringe for intraocular use, together with one sterile 23 G cannula.

**Description:** Methylvisc 2.0% is an isotonic, clear, viscoelastic and sterile solution for intraocular application. As an auxiliary device during surgery of the anterior segment of the eye, Methylvisc 2.0% maintains the depth of the anterior chamber and protects the surrounding intraocular tissue. Methylvisc contains 2.0% Hydroxypropyl Methylcellulose and supports the intraocular manipulations due to its viscoelastic and moistening properties. 1.0 mL of Methylvisc 2.0% contains 20.0 mg Hydroxypropyl Methylcellulose, Sodium chloride, Potassium chloride, Calcium chloride, Magnesium chloride, Sodium acetate, Sodium citrate, Water for injection.

Methylvisc 2.0% has a pH of 6.8 to 7.5 and osmolality of 300 to 390 mOsm/kg, similar to the aqueous humor.

Methylvisc 2.0% has a zero shear viscosity of approximately 10,000 mPa s.

Methylvisc does not cause any inflammatory or immunogenic reactions. Methylvisc 2.0% is free from preservatives. Due to its solubility in water, Methylvisc 2.0% can be removed at the end of surgery simply by irrigation and aspiration.

**Indications:** Methylvisc 2.0% serves as a volume substitute for aqueous humor during intraocular operations, such as lens extraction and insertion of an intraocular lens.

Methylvisc 2.0% maintains the depth of the anterior chamber during surgery and reduces the danger of traumatization of the corneal endothelium, the iris and the ciliary body by direct contact with surgical instruments.

**Contraindications:** Care should be taken in patients with hypersensitivity to any components in this solution.

**Directions for Use:** Methylvisc 2.0% is carefully injected into the anterior chamber using the supplied sterile single-use cannula. It is recommended for intraocular lens implantation to cover the implant and instruments with Methylvisc 2.0% immediately before surgery. This measure is to additionally protect the endothelium and the surrounding tissue. The Methylvisc 2.0% volume to be injected varies and depends on the type of surgery. As a substitute for viscoelastic lost by flow or irrigation, Methylvisc 2.0% can be injected several times. At the end of surgery, Methylvisc 2.0% must be removed completely using a suitable irrigation/aspiration device.

**Warnings and Precautions:** Residual viscoelastic not removed at the end of the surgical procedure is flushed out naturally via the trabecular system and Schlemm’s canal. There is, however, a risk of blockage of the drainage channels, which can lead to an increase in intraocular pressure. To prevent the risk of IOP, irrigation/aspiration is used at the end of the surgical procedure to remove residual viscoelastic. In case of IOP increase, treat with the preferred therapy. Patients who have been diagnosed with wide-angle glaucoma, severe myopia, diabetic retinopathy or uveitis are at particular risk of increased intraocular therapy.

- Do not reuse syringe and cannula. Any repeat usage of the syringe and cannula carries a risk of contamination and infection of the patient.
- Do not re-sterilize the pre-filled syringe and cannula.
- Do not use if package is damaged or opened.
- Do not use after the expiry date printed on the pack.
- Do not use if the syringe and cannula in accordance with accepted medical practice and applicable national, local and institutional requirements. There is no evidence concerning the safety of Methylvisc 2.0% in human pregnancy and lactation. Administration during pregnancy and lactation is at the discretion of the ophthalmic surgeon.

**Storage:** Store between 2 °C – 25 °C. Protect from light and shocks. Do not freeze.

Instrukcije za upotrebu	EN	ES	FR	HR	IT	KO
Methylvisc 2.0%	تعليمات الاستخدام	Metilvisc 2,0%	Metilvisc 2,0%	Metilvisc 2,0%	Metilvisc 2,0%	메틸비스크(Methylvisc) 2.0%

**Content:** Methylvisc 2.0% is a viscoelastic Hydroxypropyl Methylcellulose solution.

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